

# SENSI

## FootActive

Australian Premium Insoles

### Full-length Orthotic Insoles

Walking comfort for people with:



DIABETES



ARTHRITIS



SENSITIVE FEET



Plastazote cushioning layer

Flexible base with gel padding

Foot arch and metatarsal support



### Walking comfort for people with Diabetes, Arthritis or sensitive feet

Diabetic foot conditions develop from a combination of causes including poor circulation and neuropathy. Diabetic neuropathy can cause insensitivity; diminished ability to feel pain, heat and cold. Diabetics can develop minor cuts, scrapes, blisters or pressure sores of which they may not be aware, due to the insensitivity. These minor injuries may lead to ulceration and possibly even amputation!

FootActive Sensi has been designed to greatly reduce pressure and friction on the feet. The orthotic provides both support and cushioning along the entire length of the foot, thereby evenly distributing body weight and protecting sensitive areas from excess pressure and friction. In turn, this will help reduce the chance of developing foot injuries and ulceration.



FootActive Sensi may also help alleviate:

- Plantar Fasciitis • Heel Pain • Heel Spurs
- Metatarsalgia • Ball of foot pain
- Achilles Tendonitis

### Fitting & Care Instructions:

Please remove any existing insoles or lining from your shoes before inserting FootActive. If required, FootActive insoles can be trimmed to size with scissors at the toe end. You may clean the insoles with a warm soapy cloth. Do not soak them or put them in the washing machine or dryer.

Full-length insole for all types of lace-up shoes including work boots and sports shoes

